microlife[®]

ACCURATE

Comfortable pain-free
measurement

User Manual

Series 200
Pulse Oximeter

Important Product and Safety Information

Follow Instructions for Use. This document provides important product operation and safety information regarding this Pulse oximeter. Please read this document thoroughly before using the device and keep for future reference.

This non-invasive, fingertip, pulse oximeter utilizes red and infrared light with algorithm and circuitry to provide a spot check measurement of the user's pulse and oxygen saturation. It is intended for personal use in a non-professional setting by adults; for sports, training, and recreational type activities where such information IS NOT used for the monitoring/diagnosis/treatment of a medical condition.

This pulse oximeter IS intended to be used:

- For Sports, training, and recreational type activities where spot checked pulse and/or oxygen saturation information is desired.

This pulse oximeter IS NOT intended for use:

- As a medical device and/or a product supplying a user's physiological measurements related to the monitoring, diagnosis and/or treatment of a medical condition
- For continued monitoring
- With children 12 years old and under.

WARNING – Potentially hazardous situation that if not avoided may result in serious injury or death.

- This is not a medical device. It is not intended for the use in monitoring, diagnosing, or treating a medical condition (such as aiding in determining if a drug like oxygen should be administered). Likewise, this device is not for use:
 - In an electrosurgical unit (SU).
 - In an MRT or CT environment.
 - In an explosive or inflammable atmosphere.
 - In situations where alarms are required as this device has no alarms.
 - As a continuous monitor.
 - As part of a system for assessing clinical signs and symptoms.
 - With sterilizing methods; the device is not intended for sterilization.
 - During patient transport outside a healthcare facility.
- Ensure that children are supervised around this device; some parts are small enough to be swallowed and may become a choking hazard, and any attached lanyard may be of length to cause strangulation.
- Keep batteries away from children; harmful if swallowed.

- Performance of this pulse oximeter may be degraded if used outside of its specified operation temperature and humidity rating, or if stored outside of its specified storage temperature and humidity rating. Avoid storage in direct sunlight.
- Performance of this pulse oximeter may be degraded if the device is exposed to mechanical shock.
- Performance of this pulse oximeter may be degraded if the sensor lens is dirty or damaged; ensure the sensor lens is clean before using.

DO NOT:

- Use this device if you think it is damaged or anything appears unusual (e.g. sporadic operation, open packaging upon purchase).
- Dispose of the battery or device in fire; the battery may explode.
- Use a metal tool to replace the battery; doing so may short circuit the terminals.
- Attempt to recharge the battery as they are not rechargeable.
- Open/modify this device; inaccuracy and/or harmful injury may result
- Service the device when being used or when power is supplied (other than changing batteries). No other parts/components are accessible.

CAUTION – Potentially hazardous situation that if not avoided may result in minor/moderate injury, property damage, and/or damage to the device

- This device contains sensitive electronics components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g., mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy. Keep the device a minimum of 3.3 m from such devices when using or move the device to another location if interference is determined.
- When not using the pulse oximeter for extended lengths of time, remove the batteries to avoid potential battery leakage and damage to the monitor.
- Only open the pulse oximeter battery cover to replace the battery. When replacing the batteries, ensure all the batteries are replaced simultaneously to avoid battery damage and potential damage to the monitor. Microlife USA does not recommend using rechargeable batteries.
- Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations as domestic waste.

DO NOT:

- Use this device any longer than ½ hour at any given location on the body.
- Drop this device or expose it to strong vibrations; sensitive components may be affected resulting in inaccuracies and/or operational issues.

- Use the pulse oximeter outside of its specified operation temperature and humidity rating, or if stored outside of its specified storage temperature and humidity rating. Avoid storage in direct sunlight.
- Immerse in water or use abrasive agents when cleaning.

STANDARDS

In addition to the standards stated in the Instruction Manual:

- This device is compliant with the following standards:
AAMI/ANSI ES60601-1, ES60601-1-2, AAMI/ANSI/ISO 14971, ISO 80601-2-61, and
- Electromagnetic standards IEC 60601-1-2 along with FCC Part 15

Test Protocol (for laboratory use)

A functional tester cannot be used to assess the accuracy of a pulse oximeter monitor sensor as clinical testing is used to establish the accuracy of the device. However, a functional tester (Index2 FLUKE simulator version 2.1.3) may be used to measure how accurately the pulse oximeter is reproducing a specified test value

TYPE

IP22: Protected against solid foreign objects of 12.5 mm diameter and greater.
Protected against the ingress of angled falling water drops with the device positioned on a turntable. Keep Dry. 

Type BF Applied Part

SERVICE

When the device is not in use, battery/access (if applicable) may be replaced by the user. No other parts/components are accessible.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

This device does not require routine calibration or maintenance other than battery replacement.

Expected life of pulse oximeter: 5 Years

Federal Communications Commission (FCC) Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Changes or modifications to the product are not approved by

Microlife USA and could void the user's authority to operate the equipment under FCC jurisdiction.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment on and off, the user is encouraged to try to correct the interference by one or more of the following measures: 1) Reorient or relocate the receiving antenna, 2) Increase the separation between the equipment and receiver, 3) Connect the equipment into an outlet on a circuit different from that to which the receiver is connected, 4) Consult the dealer or an experienced radio/TV technician for help.

Pulse Oximeter - 200 Series

Instruction Manual

Table of Contents

1. Introduction

- 1.1. Your pulse oximeter
- 1.2. How does a pulse oximeter work?

2. Product overview

- 2.1. Product diagram
- 2.2. Inserting the batteries
- 2.3. Using the lanyard

3. How to use

- 3.1. What is blood oxygen saturation (SpO₂)?
- 3.2. Measuring blood oxygen and pulse
- 3.3. Low oxygen saturation & pulse alert

4. Error messages

5. Care and maintenance

6. Limited warranty

7. Technical specifications

8. How to contact us

1. Introduction

1.1. Your pulse oximeter

Thank you for purchasing the Microlife pulse oximeter. Your pulse oximeter is designed to provide fast and reliable non-invasive digital spot checks of your oxygen saturation of arterial hemoglobin (SpO₂) and your pulse rate. It offers clinically proven accuracy and has been designed to be user friendly.

Before using your pulse oximeter, please read this instruction manual carefully to ensure correct use. If you have additional questions regarding your oxygen saturation of arterial hemoglobin (SpO₂) and/or your pulse rate measurements, please contact your doctor.

1.2. How does a pulse oximeter work?


Fingertip pulse oximeters works by applying a mathematical formula based on Lambert Beer's Law, to measure spectrum absorption characteristics of deoxygenated hemoglobin (Hb) and oxyhemoglobin (HbO₂) in red and near-infrared zones.

This fingertip pulse oximeter works by using photoelectric oxyhemoglobin inspection technology, in accordance with capacity pulse scanning and recording technology, so that two beams of different wavelengths act as sensors (660nm red and 905nm near infrared light) when focused into a human nail tip through a finger sensor. This signal results in a measurement, which is then displayed on the screen.

2. Product overview

2.1. Product Diagram

1. ON/OFF button
2. Oxygen saturation (value as percentage)
3. Pulse rate (value in beats per minute)
4. Pulse bar
5. Low battery indicator


2.2. Inserting the batteries

After you have unpacked your device, insert the batteries. The battery compartment is located on the bottom of the device. Remove the battery cover by sliding it in the direction shown. Insert the batteries in the direction shown by the polarity symbols.

Replace the batteries when the low power indicator appears on the display.

Note:

- Use 2 x AAA 1.5 V batteries
- Always replace both batteries at the same time.

Inserting the battery


low power
indicator

2.3. Using the lanyard

To use the lanyard:

- Thread the thinner end of the lanyard through the hanging hole at the rear end of the device.
- Thread the thicker end of the lanyard through the threaded end before pulling it tightly.

Using the lanyard

3. How to use

3.1. What is oxygen saturation (SpO₂)?

Oxygen saturation indicates the percentage of hemoglobin in arterial blood that is loaded with oxygen. This is a very important parameter for the respiratory circulation system. Many respiratory diseases can result in lower oxygen saturation within human blood.

The following factors can reduce oxygen saturation:

- Automatic regulation of organ dysfunction caused by anesthesia
- Intensive postoperative trauma
- Injuries caused by some medical examinations


These situations may result in light-headedness, asthenia and vomiting. Therefore, it is very important to know the oxygen saturation of a patient so that doctors can detect problems in a timely manner.

3.2. To measuring blood oxygen and pulse


To measure your blood oxygen and pulse;

1. Insert the batteries as described in the batteries section
2. Insert one finger (nail side up; index or middle finger is recommended) into the finger opening of the device. Be sure to fully insert the finger so that the sensors are completely covered by the finger.
3. Release the device allowing it to clamp down on the finger.
4. Press the ON button

Your measurement values will appear on the screen after a few seconds, remove your finger from the device once your measurement is complete. The device will automatically switch off after approx. 8 seconds after the finger is removed from the device.

Inserting a finger

On/Off (Power) button

Notes:

- Do not shake your finger during the test. It is recommended that you stay still and do not move while taking a reading.
- The height of the pulse bar is an indication of the pulse and signal strength. The bar should be greater than 30% for a proper reading.
- The device must be able to measure the pulse properly to obtain an accurate SpO₂ measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO₂ measurement.
- Do not use the device on a single location for more than 30 minutes in order to ensure correct sensor alignment and skin integrity.

Pulse bar

Inaccurate measurements may occur if:

- The patient suffers from significant levels of dysfunctional hemoglobin (such as carboxyhemoglobin or methemoglobin).
- Intravascular dyes such as indocyanine green or methylene blue have been injected into the patient.
- Used in the presence of high ambient light (e.g. direct sunlight). Shield the sensor area with a surgical towel if necessary.
- There is excessive patient movement.
- The patient experiences venous pulsations.
- The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia.
- The patient is in cardiac arrest or is in shock.
- Fingernail polish or false fingernails are applied.

3.3. Low oxygen saturation & pulse alert

A warning alert will sound if the device detects:

- Your pulse rate is lower than 50 bpm
- Your pulse rate higher than 130 bpm
- If the SpO₂ level is lower than 94%

4. Error messages

Description	Possible Cause	Solution
SpO ₂ or pulse rate do not display normally.	1. Finger is not inserted correctly. 2. Patient SpO₂ value is too low to be measured. 3. There is excessive illumination.	1. Reinsert your finger 2. Measure 3 more times. If you determine the device is working correctly, consult your doctor.
SpO ₂ or pulse rate is unstable.	1. Finger might not be inserted deeply enough. 2. Excessive patient movement.	1. Reinsert your finger 2. Sit calmly and take a new measurement.
The device can not be powered on.	1. No batteries or low battery power. 2. Batteries are not installed correctly. 3. The device may be damaged.	1. Replace the batteries. 2. Remove and reinstall the batteries. 3. Contact customer service.
The display suddenly switches off.	1. The device automatically powers off after 8 seconds when no signal is detected. 2. The battery power is too low to operate.	1. Functioning normally. 2. Replace the batteries.

5. Care and maintenance

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the silicone that touches the finger inside of the device. Also clean the finger being tested using alcohol before and after each test. Allow the device to dry thoroughly before use.

Note:

- Never use abrasive cleaning agents, thinners or benzene (gas) for cleaning
- Never immerse the device in water or other cleaning liquids.
- Never open the device. This voids the manufacturer's warranty.

6. Limited warranty

Your pulse oximeter is warrantied for 1 years by Microlife USA Inc, against manufacturer defects for the original purchaser only, from date of purchase.

This warranty applies to the pulse oximeter only. Batteries are not covered by this warranty.

The warranty does not apply to consequential and incidental damages, or damage caused by batteries, improper handling, and accidents. Professional use, not following the operating instructions, and alterations made to the monitor or accessory by third parties, are also not included in this warranty. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

Microlife USA Inc will investigate your concern. A monitor or accessory as defined by this warranty, determined to be out of specification, will be replaced and shipped to you at no cost. A monitor or accessory as defined by this warranty, determined to be within specification, will be returned to you with a report of findings, at no cost.

Please use the below customer service contact information to reach Microlife USA Inc. regarding any warranty concerns. We ask that you please contact us before sending any product back in order to better identify, and more quickly process, your concern.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

7. Technical Specifications

Type: Fingertip Pulse Oximeter OXY 200

Display: LED display

SpO₂

Measurement Range: 70 ~ 100 %

Accuracy: 70 ~ 100 %: ± 2 %

Resolution: 1 %

Pulse rate:

Measurement range: 30 ~ 250 bpm

Accuracy: 30 ~ 99 bpm: ± 2 bpm

Resolution: 1 bpm

Operating conditions:

5 - 40 °C / 41 - 104 °F

15 - 80% relative maximum humidity

Storage conditions:

-10 - +50 °C / 14 - 122 °F

10-93% relative maximum humidity

Expected service life: 5 years

Automatic Switch-off:

Automatic shut -off in 8 seconds, when no or low signal is detected.

Battery: 2 x AAA, 1.5 V

Dimensions: 62 x 37 x 32 mm

Weight: 50 g with batteries

Technical modifications are reserved.

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